




Designed and developed to meet the requirement of ICSI procedure

Eliminating vibrations and ensuring a fast dampening effect is essential when working with inverted microscope

- Model: MS powder coated AVT/SS AVT
- The Anti-Vibration handles from light to heavy inverted microscope
- MS powder coated AVT top with stainless steel platform / In SS AVT top is frosted toughened glass platform
- The damping systems use column-friction for optimal result
- The damping systems use column-friction for optimal result
- An equilibrated and isolated platform for the inverted microscope
- Working range 15-75kg
- Roomy workspace around the micro manipulating systems
- Space for mini incubator
- MS powder coated shelf for accessories / SS shelf for accessories.
- Dimension:
1200 L x 700 W x 705 H in mm



A Product From Your Friend

Esquire IVF Equipments Pvt Ltd

#5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069
www.esquires.in | admin@esquires.in