

# HEATING BLOCK



- Made of GI with powder coat finish
- Ormon PID controller
- Temperature accuracy  $\pm 0.2^{\circ}\text{C}$
- Temperature range ambient to  $49.9^{\circ}\text{C}$
- 3 No's of anodised aluminium blocks
- Eachg blocks can hold 6 No's of ASP tubes
- Optimal  $30^{\circ}$  angle block
- Easy to clean block

*A Product From Your Friend*

Esquire Equipments

#3/16, Thiruveethi amman koil street, Korattur post, Mannurpet, Chennai 600080  
www.esquires.in | admin@esquires.in