



- Made of GI with powder coat finish
- Working Table S.S. 304
- Working area : 650 L x 600 W x 600 H in mm .
- Vertical down flow laminar system
- HEPA filters with rating of 99.99% efficacy 0.3 $\mu$  particulate
- EBM german blower for continues operation.
- Illumination light & germicidal UV Incorporated
- Additional power point inside the working area

*A Product From Your Friend*

Esquire Equipments

#3/16, Thiruveethi amman koil street, Korattur post, Mannurpet, Chennai 600080  
www.esquires.in | admin@esquires.in